



### EASY TO USE

A convenient container for solution preparation is included in the complete set. Solution for rinsing is prepared using a mixture of ordinary salt and boiled or distilled water.



### PRESSURE MONITORING

Choose the most acceptable jet pressure for you. 3 pressure levels and corresponding light indication.



### USB CHARGING

For device charging mains adapter with USB connector can be used. Wire socket on the device has a high degree of moisture protection.



### CONTAINER FOR STORING

It is convenient to store and transport the device with all accessories in a container included in the complete set.

LITTLE DOCTOR INTERNATIONAL (S) PTE. LTD.  
Yishun Central P.O. Box 9293 Singapore 917699, Fax: 65-62342197  
E-mail: info@littledoctor.sg  
www.AQUARIN.sg

© Registered trade mark of Little Doctor International (S) Pte. Ltd.  
© Copyright Little Doctor International (S) Pte. Ltd., 2023

E730-2306-1\_ENG

# AQUARIN®

OTORHINOLARYNGOLOGICAL DEVICE FOR NASAL RINSING

## NASAL RINSING AND IRRIGATION

- SAFE
- PAINLESS
- CONTROLLED




www.AQUARIN.sg


# AQUARIN®

The Aquarin device is designed for safe, controlled and painless nasal rinsing and irrigation, in order to prevent and treat inflammatory diseases of the nasal cavity in adults and children. The use of the device is recommended for prevention and treatment of inflammatory diseases of the nasal cavity at home and in prevention and treatment facilities in order to prepare patients for surgery (removal of polyps or adenoids, correction of defects and deformations of the bone and cartilage wall, etc.) and rehabilitation.

Usage of nasal irrigator helps relieve the symptoms of a cold or flu, recover from nasal congestion, moisturize the nasal mucosa and improve breathing. Nasal irrigator is especially effective in the presence of allergic reactions to flowering plants or animals.



Rinsing the openings of paranasal sinuses and nasopharyngeal cavity has a beneficial effect on the work of the respiratory system as a whole.




## FOR ADULTS AND CHILDREN

Baby nozzle is designed for familiarization and habituation of nasal irrigator usage by children. After a child gets used to the procedure and masters the skills of working with the device, it is possible to use an adult nozzle.

## COMPLETE SET

The set of the device has everything you need:

1. Nasal irrigator LD-N1
2. Adult nozzle LD-SN01
3. Nozzle for children LD-SN02
4. Tube with holder
5. Container
6. USB cable for charger
7. Instruction Manual
8. Warranty card



## CHARACTERISTICS

Model	LD-N1
Rated voltage/current	5V / 500mA
AB capacity	2x1300 mA/h
Water pressure: position «●»	minimum
position «●●»	middle
position «●●●»	maximum
Service life* - irrigator	5 years
- AB	about 3 years
- nozzle	about 1000 cycles
The recommended time of continuous use is	no more than 5 minutes
The time required for the device to cool	is at least 10 minutes
Container capacity	not less than 900 ml
Liquid consumption, approximately, ml/min	100 - 500, depending on the position of the MODE switch
Enclosure protection degree	IPX7
Overall dimensions of the device, mm	(154±2)×(58±2)×(51±2)
Weight (without packaging), g	470±10